# Antimicrobial Susceptibility Test System Delivering Phenotypic MICs in Hours from Positive Blood Cultures for Fastidious- and Non-Fastidious Pathogens


## Introduction

Rapid identification (ID) and antimicrobial susceptibility testing (AST) direct from positive blood cultures (BCF) are important in helping provide septic patients with earlier and better optimized antimicrobial treatment. The great variety of different pathogens isolated from positive BCFs puts high demands on the identification technology, and around 10% of all pathogens isolated from BCFs are so-called fastidious bacteria (require growth supplements). Currently, no 'direct AST from positive blood culture' system for fastidious bacteria is available.

We have evaluated the ability to perform AST directly on a positive BCF for both fastidious and non-fastidious pathogens and report the MIC after just six-hour's growth. Our method comprised automated sample preparation and image-based antimicrobial susceptibility testing in both standard CAMHB as well as CAMHB with additives to support the growth of fastidious pathogens.

## Materials and methods


The AST results presented are obtained using a prototype system. Blood culture flasks spiked with blood from healthy individuals and clinical isolates were cultured in a Bactec 9050 cabinet until signaled positive. A 500 µl aliquot from each positive flask was subjected to automated sample preparation followed by inoculation into discrete concentrations of different antimicrobials at several dilution ranges. Time-lapse images were taken over 18 hours. MIC is determined as the lowest concentration of antibiotic inhibiting growth (Fig 1). The MIC determined by a proprietary automated algorithm from time-lapse images after six hours growth was compared to the MIC determined manually via visual observation after 18 hours growth. If the same isolate was used in several independent spike experiments, the median value of the 18-hour MICs from that strain was used as reference MIC.



**Fig 1.** An aliquot from a positive blood culture is subjected to automated sample preparation and adjustment, transferred to a consumable with dried antibiotics for broth microdilution AST, and imaged with time-lapse imaging. The MIC is calculated from the images based on the lowest concentration of antibiotics inhibiting growth measured as total biomass in the image.

#### www.qlinea.com

Q-linea and ASTrID are registered trademarks and ASTar is a trademark of Q-linea AB. Products are for research use only. © Q-linea AB 2018. Q-linea AB, Dag Hammarskjölds väg 52 A, SE 752 37 Uppsala, Sweden. Markus Klintstedt, Ylva Molin, Magnus Sandow, Lotta Levén, Jan Grawé, Charlotta Göransson, Jenny Göransson, Mats Gullberg. Q-linea, Uppsala, Sweden.



## Results

Both fastidious and non-fastidious bacteria can be analyzed using time-lapse microscopy (Fig 2). The MIC values obtained after six hours growth achieved an Overall Essential Agreement of more than 95% compared with MICs determined after 18-hour's growth (Table 1 and 2). Approximately 30% of the different bacteria-antibiotic combinations tested had a resistant or intermediate phenotype. To challenge the system, around 20% of all bacteria-antibiotic combinations had an MIC for that antibiotic on a clinical breakpoint. This explains why the Overall Categorical Agreement is lower than the Overall Essential Agreement obtained in the overall data set.




Fig 2. Biomass plotted over time for a selected number of bacteria and antibiotics. A growth / no growth assessment is made per concentration of each antibiotic to determine the MIC.

**Table 1.** Summary of Overall Essential Agreement (EA) and Overall Categorical Agreement (CA) for six-hour image-based AST data of fastidious pathogens isolated from blood cultures spiked with blood from healthy individuals and clinical isolates and grown to positivity. (total N = 124) mE = minor error, ME = major error, VME = very major error.

Streptococcus pneuomniae	Total no. of tests	EA (%)		No. of tests						
			<b>CA (%)</b>	S		R	mE	ME	VME	
Ampicillin	20	100%	80%	6		14	4			
Vancomycin	17	100%	100%	17						
Cefotaxime	19	100%	89%	6	7	6	2			
Streptococcus mitis										
Ampicillin	11	90%	100%	11						
Vancomycin	12	100%	100%	12						
Benzylpenicillin	12	100%	73%	6	6		2			
Haemophilus influenzae										
Cefotaxime	11	91%	91%	1		10			1	
Ciprofloxacin	11	100%	100%	8		3				
Meropenem	11	82%	100%	11						

**Table 2.** Summary of Overall Essential Agreement (EA) and Overall Categorical Agreement (CA) for six-hour image-based AST data of pathogens isolated from blood cultures spiked with blood from healthy individuals and clinical isolates and grown to positivity (total N = 672) mE = minor error, ME = major error, VME = very major error.

	Total no. of tests	EA (%)			No. of tests					
Antimicrobial agent			CA (%)	S		R	mE	ME	VME	
Amoxicillin-clavulanic acid	27	96%	100%	9		18				
Ampicillin	30	97%	100%	27		3				
Benzylpenicillin	6	100%	100%			6				
Cefotaxime	48	100%	81%	22	3	23	9			
Cefoxitin	10	90%	80%	6		4			2	
Ceftazidime	66	95%	86%	43	6	17	9			
Ceftolozane-tazobactam	16	75%	100%	12		4				
Ciprofloxacin	81	99%	96%	43	3	35	3			
Colistin	38	92%	100%	26		12				
Daptomycin	11	100%	100%	8		3				
Gentamicin	64	98%	95%	39	3	22	3			
Levofloxacin	14	100%	100%	3		11				
Meropenem	81	91%	100%	78	3					
Piperacillin-tazobactam	67	94%	84%	61	6		6	5		
Tetracycline	14	100%	93%	4	3	7	1			
Tobramycin	41	93%	88%	35		6	5			
Frimethoprim/sulfamethoxazole		100%	100%	8		3				
Vancomycin	47	100%	94%	41		6			3	
Fotal	672	95.7%	93.1%	465	27	180	36	5	5	

Gram negative: P. aeruginosa, E. coli, A. baumanni, P. mirabilis, K. pneumoniae, E. cloacae Gram positive: S. aureus, E. faecalis, S. epidermidis

### Conclusions

This study shows that it is possible to determine MICs for fastidious pathogens isolated directly from a positive BCF using image-based AST. Time-lapse images show that a common algorithm can detect differential growth within six hours for a large variety of bacteria / antibiotic combinations. MIC was determined by extracting biomass information in each concentration of antibiotic. The imaging is compatible with the additives that support the growth of fastidious organisms and the variable effects on bacterial morphology caused by a broad range of different antibiotics.

Overall EA: 95.7% CA: 93.1%